

CLEVER SOLUTIONS AT AFFORDABLE PRICES



GynoPrep[®]
THIN-LAYER PREPARATION
FOR *Cellspin*[®] CYTOCENTRIFUGES

> Compatible with Shandon™ Cytospin™



GynoPrep®

The successful medium in the Cervical cancer screening


With the appearance of gynecological cancer prevention, the methods developed continually. Better technology means higher precision, which in turn means better safety for the patients.

GynoPrep® contributed as one of the leading diagnostic procedure, that cervical cancer become an exception today.

GynoPrep® meets the chemical requirements as a carrier medium for HPV typing according to the new guideline for gynaecological cancer screening. In addition, it also meets the requirements of the European Directive 98/78/EC and is CE IVD certified.

GynoPrep® by comparison with standard Pap smears:

Standard Pap smear




The **standard Pap smear** shows disadvantages:

- › Not all cells were investigated
- › Cells can be overlap
- › Blood can overlap the cells

All this complicates the assessment.

Smear with GynoPrep®



The smear with **GynoPrep®** is an alternative to the normal smear:

- › Clean and clear smears
- › Preparations can be better diagnosed by the doctor
- › Smear will be extracted with a special sampling brush. Extraction is absolutely painless
- › After extraction the brush will be put into the **GynoPrep®** liquid then the liquid will be given to the laboratory
- › So approx. 100 % of the cells will arrive in the laboratory
- › It allows a more exact assessment

All this is especially important to detect and treat preparatory stages of cervical cancer at an early stage.

GynoPrep® can be processed manually and also runs in Cobas automatic machines. Choose our kits from disposable or reusable systems.

A **GynoPrep®**-Kit contains all necessary components - from the collection to the finished preparation.

GynoPrep® can also be ordered in a single carrier vessel.

GynoPrep® is also used by the German Armed Forces.

GynoPrep® - Kit

Item no.	Item
330-B	GynoPrep®-Kit - Rovers® Cervex-Brush® Combi 100x Rovers® Cervex-Brush® Combi Brushes 100x GynoPrep® fixative 100x ECOseal 100x ECOslides (slides), coated 2x GynoPrep® dilution fluid, 30 ml
330-B-D (ECOfunnel®, disp.)	GynoPrep®-Kit - Rovers® Cervex-Brush® Combi 100x Rovers® Cervex-Brush® Combi Brushes 100x GynoPrep® fixative 100x ECOfunnel®, disposable 100x ECOslides (slides), coated 2x GynoPrep® dilution fluid, 30 ml
330-S (individ. wrapped)	GynoPrep®-Kit - Rovers® Cervex-Brush® Combi (individually wrapped) 100x bags for single use with: 1x Rovers® Cervex-Brush® Combi Brush 1x GynoPrep® fixative 100x ECOseal 100x ECOslides (slides), coated 2x GynoPrep® dilution fluid, 30 ml
450-B	GynoPrep®-Kit - Rovers® Cervex-Brush® Combi 100x Rovers® Cervex-Brush® Combi Brushes 100x GynoPrep® fixative 100x ECO filter cards 100x ECOslides (slides), coated 2x GynoPrep® dilution fluid, 30 ml

GynoPrep® - Vials

Item no.	Item
351-2	GynoPrep® - 100x vials, filled with 20 ml fixative

GynoPrep® - Accessories

Item no.	Item
332-100	GynoPrep® - Rovers® Cervex-Brush® Combi brushes - PU 100 pcs.



THARMAC® GmbH Laboratory Solutions

Commercial unit

Hasselborner Straße 19-21 » DE-35647 Waldsolms » PHONE: +49 (0) 6085 98 99 1-0
FAX: +49 (0) 6085 98 99 1-19 » MAIL: info@tharmac.de

Place of business

Borsigstraße 7A » DE-65205 Wiesbaden » PHONE: +49 (0) 6122 5 88 97-33
FAX: +49 (0) 6122 5 88 97-36 » MAIL: info@tharmac.de